

Microbiological Drinking Water Analysis

Test	Method	Bottle	Perservative	Hold Time	Turn Time
Total Coliform MMO-MUG	SM 9223-18	100ml Sterile	Sodium Thiosulfate	30 hours	36 hours
E-Coli MMO-MUG	SM9223 B	100ml Sterile	Sodium Thiosulfate	30 hours	36 hours
Total Coliform Quanti Tray	SM9223 B-QT	100ml Sterile	Sodium Thiosulfate	30 hours	36 hours
Total Coliform Membrane filter	SM 9222 G	100ml Sterile	Sodium Thiosulfate	30 hours	96 hours
Fecal Coliform Memebrane filter	SM 9222 D	100ml Sterile	Sodium Thiosulfate	30 hours	36 hours
E.Coli Membrane Filter	SM 9223 D	100ml Sterile	Sodium Thiosulfate	30 hours	96 hours
Total Plate Count	SM 9215 B	100ml Sterile	Sodium Thiosulfate	30 hours	96 hours